EURODENT

Treat yourself!

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Intended use of the device

The ISOFLEX dental unit is a medical device intended for performing dental treatments with the aid of specially installed accessories.

Environment of use

The ISOFLEX dental unit is intended for use in health care environments, such as hospitals, outpatient clinics or dental practices.

Qualified personnel

The device and its accessories may only be used by professional medical personnel specialised in dentistry who:

- are authorised to use the device
- are aware of the risks involved in using the equipment
- are aware of the standard procedures for performing the treatments for which the ISOFLEX dental unit is intended.

Limits of use

The device may not be used for operations outside the defined intended use. The user is responsible for assessing the suitability of the physical and clinical condition of the patient to be treated.

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Technical data

Trade mark: EURODENT
Model: ISOFLEX
Manufacturer: PROMED Srl
 Via Rossini, 26
 40050 Villanova di Castenaso
 (Bologna), Italy
 Tel: +39 051 0390446
 Fax: +39 051 781 625
 E-mail: info@promeditaly.it
Rated voltage: 230 V
 (~ alternated)
Rated frequency: 50 Hz
Rated power: 460 VA
Weight: 130 Kg
Maximum lifting load of the chair: 280 Kg
Classification: Class I based on the electrical protection level; Type B based on the level of direct or indirect electrical insulation.
 Classification items based on the EN 60 601-1 standard:
 5.3 Common appliances (appliances whose casing is not protected against water penetration);
 5.5 Appliances not suitable for being used in the presence of a flammable anaesthetic mixture with air, oxygen, or nitrous oxide;
 5.6 Appliances for continuous operation with intermittent loading.

Dimensions:

The structure of the dental unit is custom-made to optimise the footprint and for the best possible usability.
 An outline diagram and the available degrees of freedom are shown below.

The equipment is conformed to EEC directive 93/42/CEE "Medical Devices"

The manufacturer accepts the responsibility for safety properties of our equipment only if:

- Installation, maintenance, repair and modifications are carried out by a person explicitly authorized by us;
- The general and electric installation of the room is conformed to the related laws;

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- The equipment has been correctly set up according to our "operating instructions";
- The equipment has been correctly used according to our "User Manual".

The manufacturer declares that the equipment is a complex dental unit for being used on the man only by a high qualified practical staff in dental practice as job for therapy.

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Warning: any modification to this equipment is prohibited

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Warning: to avoid the risk of electrical shock, this unit must be connected exclusively with a network electrical connecting provided of grounding wires circuit.

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We recommend checking the electrical safety parameters, leakage currents, insulation, and earth conductor continuity every year by an authorised technician.

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WARRANTY, CE MARKING, AND TRACEABILITY OF THE PRODUCT (EUROPEAN DIRECTIVE 93/42)

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The manufacturer ensures that the device has been built correctly, using top quality raw materials, carrying out all the tests required, and complying with all the standards in force. The warranty of this product starts from the date of delivery, which must be proven by returning the registration card, duly filled out and signed by the purchaser.

The warranty is valid under the conditions specified in the warranty booklet supplied with the product. The warranty becomes automatically null and void if the product is installed, repaired, modified, or tampered with by the purchaser or unauthorised third parties; if the

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Registration Card, duly filled out and signed, was not sent to the manufacturer immediately after the installation of the device.

Only refer to authorised support centres or the manufacturer for installing the device, receiving support, and obtaining original spare parts. The responsibility for installation, maintenance, and repair interventions is always of the support centre that has carried out these operations.

The manufacturer is not liable for damage resulting from interventions of any kind carried out by unauthorised personnel and/or the use of non-original spare parts.

In compliance with the regulation in force and in order to keep the CE Marking (Dir. 93/42/EEC), the purchaser and/or the dealer must fill out and return the product Registration Card upon the installation of the device, thereby allowing the manufacturer to ensure the product traceability procedure.

Warning! Both the warranty and CE Marking are automatically null if the Registration Card is not returned to the manufacturer duly filled out within the required terms and in case of interventions carried out by unauthorised people.

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APPLIED PARTS

Based on its class, protection rating, and electrical insulation, this device is suitable for being used with type B  or BF .

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The manufacturer will not check non-original spare parts and accessories. Therefore, their compliance, safety, and performance cannot be certified.

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ESSENTIAL PERFORMANCE

This medical device supports the patient in a position suitable for the dentist's interventions. Electromagnetic disturbance does not affect the essential performance of the device.

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SAFETY DEVICES

The chair and the dental unit are provided of many devices to ensure operating safety.

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INSTALLATION

The medical device is designed to be exclusively installed in dental practices, clinics and hospitals, i.e. in professional healthcare facilities (*Professional Healthcare Facility Environment*).

Make sure that you have enough space for the dental unit; The authorised service centres can provide you with the installation template in a 1:1 scale, which is also provided with every unit sold.

The electrical system of the facility in which the device is installed must comply with the safety standards in force concerning electrical systems in medical premises.

PROMED will not be held liable for damage to persons or objects caused by the failure to comply with this rule.

Only use a suitable power supply source and connection, as indicated on the identification plate (placed under the chair's seat). Only refer to authorised support centres or the manufacturer for installing the device, receiving support, and obtaining original spare parts.

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CHAIR WORKING CYCLE (DUTY CYCLE)

Continuous operation: 2 minutes - Pause (stand-by): 18 minutes.

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ELECTROMAGNETIC INTERFERENCE

This device complies with standards concerning electromagnetic emissions. However, despite complying with the regulations in force, interference may still occur on any installation. In this regard, we recommend switching off mobile phones and/or any type of transmitter when you are in the room where the dental unit is installed.

Do not connect electrical equipment to the dental unit, including other electrical cables not mentioned in

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the list of cables provided in this manual.

Cables and accessories not included in this list may cause electromagnetic interference. Position any equipment (e.g., electrosurgical scalpels, suction devices, etc.) as far away as possible from the electronic boards of the dental unit, or, in any case, at least 1 metre away from it.

Connect any foreign equipment to sockets on circuits other than those used for the dental unit.

Contact the manufacturer or an authorised support service, if needed.

Any anomaly resulting from interference can be detected by repeatedly switching on and off the source of the interference.

To know the electromagnetic interference test limits, refer to relevant chapter of this Instruction, which contains the results of the tests carried out in compliance with the standards in force.

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INTENDED AND UNINTENDED USE

This medical device is intended for outpatient dental treatments.

Do not open the device's casing and carry out only the adjustments shown in the operating instructions. Contact the authorised service if the dental unit is not working. The manufacturer will not be held liable in the following cases: assemblies, calibrations, changes, replacements, and repairs carried out by unauthorised personnel; use of non-original spare parts; electrical system upstream the device is not conforming to the standards in force; the device is not installed, repaired, and used in compliance with the instructions provided in this manual.

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(*) *ISM (industrial, scientific and medical) bands* between 0,15 MHz and 80 MHz are: 6,765 MHz to 6,795 MHz; 13,553 MHz to 13,567 MHz; 26,957 MHz to 27,283 MHz; and 40,66 MHz to 40,70 MHz *Amateur radio bands* between 0,15 MHz and 80 MHz are: 1,8 MHz to 2,0 MHz, 3,5 MHz to 4,0 MHz, 5,3 MHz to 5,4 MHz, 7 MHz to 7,3 MHz, 10,1 MHz to 10,15 MHz, 14 MHz to 14,2 MHz, 18,07 MHz to 18,17 MHz, 21,0 MHz to 21,4 MHz, 24,89 MHz to 24,99 MHz, 28,0 MHz to 29,7 MHz and 50,0 MHz to 54,0 MHz

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Electromagnetic Compatibility

PERFORMANCE EMC CEI EN 60601-1-2	
The ISOFLEX dental unit is designed to be used in the electromagnetic environment specified below. The client or user of the ISOFLEX dental unit must ensure that it is used in such an environment.	
 WARNING: the use of accessories, transducers and cables other than the specified items, except those sold by the manufacturer as spare parts, may cause an increase in emissions or a decrease in the immunity of the device.	
Emissions EMC – Class CISPR 11	Group 1
Emissions EMC – Harmonic Class (EN 61000-3-2)	Class A
Immunity EMC – ESD (according EN 61000-4-2)	15 kV in air; 8 kV in contact
Immunity EMC – burst (according EN 61000-4-4)	2 kV 100kHz (L-N-PE)
Immunity EMC – surge (according EN 61000-4-5)	1 kV -100 kHz (L-N); 2kV 100kHz (L-PE, N-PE)
Immunity EMC - magnetic field (according EN 61000-4-8)	30 A/m
Immunity EMC - RF fields (according EN 61000-4-3)	3 V/m from 80MHz to 2,7 GHz modulation 1kHz AM 80% 27 V/m 385MHz modulation 18Hz PM 50% 28 V/m 450MHz modulation FM 1kHz sine ±5Hz PM 50% 9 V/m 710-745-780 MHz modulation 217Hz PM 50% 28 V/m 810-870-930 MHz modulation 18Hz PM 50% 28 V/m 1720-1845-1970 MHz modulation 217Hz PM 50% 28 V/m 2450 MHz modulation 217Hz PM 50% 9 V/m 5240-5500-5785 MHz modulation 217Hz PM 50%
Immunity EMC - RF currents (according EN 61000-4-6)	3 Veff 150kHz-80MHz modulation 1kHz AM 80% 6 Veff ISM Band and Amateur Radio Band (*)

(*) ISM (industrial, scientific and medical) bands between 0,15 MHz and 80 MHz are: 6,765 MHz to 6,795 MHz; 13,553 MHz to 13,567 MHz; 26,957 MHz to 27,283 MHz; and 40,66 MHz to 40,70 MHz Amateur radio bands between 0,15 MHz and 80 MHz are: 1,8 MHz to 2,0 MHz, 3,5 MHz to 4,0 MHz, 5,3 MHz to 5,4 MHz, 7 MHz to 7,3 MHz, 10,1 MHz to 10,15 MHz, 14 MHz to 14,2 MHz, 18,07 MHz to 18,17 MHz, 21,0 MHz to 21,4 MHz, 24,89 MHz to 24,99 MHz, 28,0 MHz to 29,7 MHz and 50,0 MHz to 54,0 MHz

List of cables that can be used with this dental unit

Code	Description	Fuction	Note
D1021860	PATIENT CHAIR WIRING	Patient chair board	
D1021870	WHITE TURBINE WIRING	Turbines	
D1021880	YELLOW ULTRASONIC SCALER WIRING	Ultrasonic scaler	
D1021890	BLUE MOTOR WIRING	Motor	
D1021900	DASSYM EU-33-LFA MOTOR WIRING	Motor	
D1021910	DASSYM SERIAL WIRING	DASSYM motor	
D1021920	MCX MOTOR WIRING	STD&CP MCX MOTOR board	
D1021930	MCX SERIAL WIRING	MCX STD&CP	
D1021940	SATELEC ULTRASONIC SCALER WIRING	SATELEC ultrasonic scaler	
D1021950	WATER SPRAY SWITCH WIRING	Water spray wiring	
D1020390	TURBINE-MOTOR SOLENOID VALVE WIRING	Turbines/Motors - no CP	
D1020470	ULTRASONIC SCALER SOLENOID VALVE WIRING	Ultrasonic scaler	
D1340110	TURBINE CONNECTOR	Turbine	
D1340120	MCX MOTOR CONNECTOR	STD&CP MCX MOTOR board	
D1340130	DASSYM MOTOR CONNECTOR	DASSYM motor	
D1340140	ULTRASONIC SCALER CONNECTOR	Ultrasonic scaler	
D1021600	PEDAL CABLE	Pedal cable	
C0690004	6-PIN MICRO MATCH WIRING	Pedal joystick	
D1021840	PEDAL W/ACC REED SENSOR WIRING	Pedal wiring	
D1021850	PEDAL CONNECTOR W/ACC WIRING	Pedal wiring	
D1021960	ACCELERATOR PEDAL BUTTON WIRING	Pedal wiring	
D1021970	PEDAL CABLE W/ACCELERATOR	Pedal cable	
D1021590	PEDAL WIRING IN PATIENT CHAIR - EVO L=800	Patient chair/Pedal cable connector board	
D1020290	SWIDENT EARTH WIRING	Patient chair board	
D1021400	POTENTIOMETER WIRING	Patient chair/Limit stop potentiometers board	
D1020270	SWIDENT SWITCH WIRING	Patient chair/Patient chair switch board	
D1020280	PATIENT CHAIR SAFETY WIRING	Patient chair board	
D1020310	GENERAL WATER SV WIRING	Patient chair/EV H2O board	
C0940002	POWER WIRE SCHUKO 1.5 m 3x0.75 BLACK	Patient chair power supply	

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
(*) ISM (*industrial, scientific and medical*) bands between 0,15 MHz and 80 MHz are: 6,765 MHz to 6,795 MHz; 13,553 MHz to 13,567 MHz; 26,957 MHz to 27,283 MHz; and 40,66 MHz to 40,70 MHz Amateur radio bands between 0,15 MHz and 80 MHz are: 1,8 MHz to 2,0 MHz, 3,5 MHz to 4,0 MHz, 5,3 MHz to 5,4 MHz, 7 MHz to 7,3 MHz, 10,1 MHz to 10,15 MHz, 14 MHz to 14,2 MHz, 18,07 MHz to 18,17 MHz, 21,0 MHz to 21,4 MHz, 24,89 MHz to 24,99 MHz, 28,0 MHz to 29,7 MHz and 50,0 MHz to 54,0 MHz

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Disposal pursuant to Legislative Decree 27 of March 2014 "Transposition of Directive 2011/65/EU on the restriction of the use of certain hazardous substances in electrical and electronic equipment" and 49 of March 2014 "Transposition of Directive 2012/19/EC on waste electrical and electronic equipment (WEEE)".

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Equipment for intermittent operation (chair: 2 minutes on - 18 minutes off)

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
Transportation and storage conditions

When the equipment is packed for shipment and/or storage, it can support the following conditions:

- Environmental temperature from -40°C to +70°C;
- it should not be exposed to rain or excessive humidity.

The other symbols on packaging are:

- maximum the number of packages stackable (3)
- the high side of the package
- the fragility of the content.

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Temperatura ambientale da +10°C a +35°C.

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Pressione atmosferica: da 700 a 1600 mbar.

Operating conditions

Environmental temperature from +10°C to +35°C;

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In addition, devices are present on the front of the dental unit attached to the back of the patient chair that indicate its status:

- Signal LED button (A).
This activates the water spray from the pressurised bottle.
The bottle is under pressure when the LED is lit.

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Warning: Do not disassemble the bottle if it is under pressure.


Fig. 5 • Pct.54 • Abb. 5

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- Water spray on indicator light (B).
- Reverse rotation warning light (C).
- LED on indicator light (D).

The LEDs all lit from left to right, during the first phase, indicate that the instrument is active. Using lever L1 will indicate the speed or power of the instrument in use by illuminating the series of LEDs up to a maximum of 5 when the instrument reaches maximum intensity.

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Warning: The functions described below are only possible with the instruments at rest, except for the syringe, which may be raised. The patient chair movements are inhibited with an instrument in operation, except for the syringe.

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
Use of the unit

Unit main switch

The unit gets ready for the operating pressing the main switch (position 1). The main switch will be illuminated by a green light.

Protection fuses

The dental unit is provided of 2 protection fuses for the power supply. To replace these fuses, you must unscrew the fuse holder (n° 2) and replace the fuse (n° 1).

Use only fuses having the same characteristics of the original ones (T3,15A - 250V - 5x20).

Programming of the patient chair positions

All keys control of the patient chair (table's dentist panel, table's assistant panel and foot control if present) have a double function: if they are pressed (>1 second), they start the driven movements; if they started by an impulse, they started the program.

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Patient chair movements and program recall (lever pedal)

Using the joystick (L3)

- Position (1) ▲ the maintained contact activates the ascent of the patient chair. One pulse activates the automatic programmed movement of the patient chair stored position 1. The lights is on when the chair stops.

- Position (2) ▼ the maintained contact activates the descent of the patient chair. One pulse activates the automatic programmed movement of the patient chair stored position 2. The lights is on when the chair stops.

- Position (3) ◀ the maintained contact activates the descent of the backrest. One pulse activates the automatic programmed movement of the stored position of chair 3, which returns it to zero. The lights is off when the chair starting moving.

- Position (LP) ▶ the maintained contact activates the ascent of the backrest. One pulse activates the automatic programmed movement of the stored position of chair 4, which positions it at rinsing height with cup filling, timed cuspidor bowl washing and then return to the last working position. The light is on after the second pulse.

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Other controls with instruments at rest or syringe activated.

- Switching the light on and off. To switch the lamp on or off, turn the lever (L1) to the end and release it.
- Assistant call bell activation Fully press the pedal (L2) to activate the bell.

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Tool-activated commands

Warning: The functions described below are only possible with only one instruments active, except for the syringe, which may be raised. The patient chair movements are inhibited with an instrument in operation.

- Instrument rotation reversal. Rotation is reversed by moving the joystick (L3) to position 2
- Water spray activation/exclusion. The water spray is activated and/or excluded by moving the joystick (L3) to position 1
- Chip activation The air blow or chip is activated by fully pressing the chip lever (L2)
- Instrument full scale adjustment **+/-** The full scale of the instruments is adjusted by activating the instrument and tapping the joystick (L3) to the + position to increase it or tapping the joystick (L3) to the - position to decrease it.

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- Switching the light on and off. 
To switch the light on or off, fully press the pedal (L1) and release it.
- Assistant call bell activation 
Press the button (5) to activate the bell.

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- Instrument rotation reversal. 
The rotation is reversed by fully pressing the button (2).
- Water spray activation/exclusion. 
The water spray is switched on and/or off by fully pressing the button (1).
- Chip activation.
The air blow or chip is activated by fully pressing the button (5)
- Instrument full scale adjustment **+/-**
The full scale of the instruments is adjusted by moving the instrument and repeatedly and fully pressing button + to increase it or repeatedly and fully pressing button - to decrease it.

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Memorization of programs

For changing one of the other four of the pre-selected functions stored (program no. 1, 2, 3 and LP), bring the patient chair into the desired position manually by pressing the program keys; then press towards the high the lower cover of the arm (Pantograph) of the patient chair and in the same time press the wanted program key and hold depressed until hearing a Beep sound, this is the confirmation of the new storage position instead of the old one.

The new stored program can now be called by pressing the appropriate key.

It is not allowed to move the patient chair when the instruments are taken from holders and are activated by the multifunctional foot control for the keys commands.

DO NOT set up programs when movements are at end stop.

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Physiological solution (Optional)

If the bottle for physiological solution is present, you can select by the key board the liquid to the instrument's destination. The liquid will be dispensed from the bottle with the button light on.

Fig. 12 • Pct. 12 • Abb. 12

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The LED indicator must be switch off during the operating substitution of the bottle.

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Water spray management system.

The instrument water spray activation is managed using the pedals as described in the "USE OF THE UNIT" section

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The manufacturer supplies these instruments with original instructions for use. Follow the instructions in the documentation supplied with the devices for technical characteristics, performances, maintenance procedures and for proper use.

The pressure regulation of the instruments

Every instrument has an individual regulation of the spray; in the lower part of the dentist's instrument table are placed the regulation taps with the following functions:

- 1 - Screw tap: it regulates the volume of air or the cooling. This adjustment must be made by a qualified technician;
- 2 - Tap with Knurl: it regulates the water spray;
- 3 - Screw tap: it regulates the air of chip blower and spray. This adjustment must be made by a qualified technician.

Dentist's syringe and assistant syringe

This instrument does not have the program functions.

- 1 - Control Key for providing water;
 - 2 - Control Key for providing air;
- Pressing both keys (1 and 2) at the same time it will provide spray.

The syringe's cover can be unthread for permitting the sterilization, pressing the realising push button 3.

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Turbine

When the instrument is taken from its holder, the LED indicates the function's power percentage, automatically it will be appeared as 100% value and it cannot be changed. If the turbine has F.O., it's switched ON (and OFF, resp.) by pressing key board (LED indicator illuminated). The instrument has 3 pressure's regulators situated in the lower part of the back side of the instrument's table.

Micromotor

When the instrument is taken from its holder, the same is arranged for the operation. The speed regulation depends on the pressure exercised on the foot control and the key of selection of the spin sense will be active. If the Fibre optic is arranged it will be active the key for switch it on / off. The instrument has the regulator of the spray placed in the lower part of the dentist's instrument table.

The air motors supplied with the unit are equipped by the standard connection type "E", conformed with a directive ISO 3964. It's expected the utilisation of the hand pieces exclusively CE marked.

Electrostatic discharges on the foot control base, can create the stop of the micromotor; in this case it is possible to restore the normal function resetting to zero the foot control (bringing lever 1 of the foot control completely to the left side).

Ultrasonic Scaler

When the instrument is taken from its holder, the LED indicates the previously memorized power and it can be changed. The instrument has one pressure's regulator for cooling water in the lower part of the doctor's table.

The oscillation of the ultrasonic scaler to the high frequency can compromise (interfered) the working of the pace-maker.

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Curing Light Mini LED

- 1 - Navigation key for select the working modality
- 2 - Key on/off of working

3 - LED indicator "progressive"

4 - LED indicator "at impulses"

5 - LED indicator "fast"

6 - LED indicator working conditions:

Green - normal

Red - Problem of electric alimentation

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Menu "progressive" mode:

select with the navigation key. In "progressive" mode (polymerization similar to menu step with halogen lamp) you will have at disposition:

- Progressive start for 10 seconds ones;
- Full power for 10 seconds one.

Menu "to impulses" mode:

select with the navigation key. This menu corresponds to the full power of the lamp in the mode to impulses, that is an emission of 10 successive flashes with a pause period of 250 Ms in the luminous exposures.

Menu "fast" mode:

after removing the instrument from its holders, the lamp is situated in fast mode. This menu activates the maximum luminous power (according to the terminal part) for 10 seconds.

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For the correct utilisation of the Curing light, it's important to act in the following way:

1. To choose the menu and the lamp is ready for using;
2. Put the terminal part nearest to the surface of the polymerization material to polymerization. It does not need that the terminal part should be in contact with the composite material, as it can be damaged considerably and to decreasing the effect (for the dirt);

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3. With one short pressing the key *ON/OFF*, it is switched the polymerization cycle. This action will be confirmed by a "beep" sound.
4. A simple pressure of the key *On/Off* is however possible to interrupt the cycle in any time.
5. If it is not used for 3 consecutive seconds, the curing lamp goes in stand-by and all indicators switched off. It is enough to sufficient press one key for switching it on.

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Optical details

Wavelength: 420-480 nm
 Intensity: from 1250 to 2000 mW/cm²
 ±10% (it depends from the final part)

Tension breaking in the power supply can cause the stopping of the polymerized lamp; in this case it is possible to restore the normal use resetting to zero the foot control (hanging the lever 1 of the foot control to the left completely).

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For the application of others hand pieces, after the first unit's supply, please contact Eurodent.

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- Lampada alogena EDI
- Lampada a led MAIA
- Lampada a led ALYA
- Lampada a led ALYA THEIATECH
- Lampada a led EVA
- Lampada a led EVA CAM
- Lampada a led EVA THEIATECH

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Dental light

The following Faro dental light models can be installed on the unit:

- EDI halogen light
- MAIA LED light
- ALYA LED light
- ALYA THEIATECH LED light
- EVA light
- EVA CAM light
- EVA THEIATECH light

Refer to the original use and maintenance manuals supplied with the devices.

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Never use the curing light without the special protection shield on and pay attention not to direct the light beam into eyes.

Optical Fibre

To insert the optical fibre onto the handpiece, push it fully down into its housing until a click can be heard.

The optical fibre is made in glass and is therefore fragile: in case of knocks or falls it may break or chip, producing sharp bits which may cause injuries and compromising final quantity of emitted light.

The lamp is equipped with permanent memory so that at power-up the last used cycle is recalled.

The lamp is also equipped with a buzzer which emits a beep when light is turned on, every 5 seconds of operation, and when light is turned off.

Furthermore, the lamp is equipped with a thermal protection.

Sequence of operations

Extract lamp from its holder on the dental unit.

Rotate lamp and/or optical fibre to the best setup for use (straight, pistol or intermediate position).

Select the cycle to be used by means of MODE button.

Put optical fibre into the appropriate position: the output of the optical fibre must be kept as close as possible to the material to be cured, though without touching it.

Start the cycle by pressing the START button and wait until light emission turns off automatically.

MODE button

Press this button to select the cycle to be used. The set cycle is displayed by the turnon of the corresponding LED (or pair of LEDs). Cycle selection is only possible when lamp is not emitting light.

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START button

Press this button to start the set cycle. If this button is pressed again during the cycle, light emission will be immediately stopped.

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FAULT SIGNALS

In case of malfunction of the curing lamp, the following signals are shown on the keyboard:

- LED 5 and LED 1 on: Light beam not being emitted by the lamp. Contact Technical Service.
- LED 5 and LED 2 on: Microcontroller for starting the lamp is malfunctioning. Contact Technical Service.
- LED 5 and LED 3 on: Insufficient power supply. Check supply voltage.
- LED 5 and LED 5 blinking: Handpiece thermal shutdown has tripped. These LEDs will continue to blink until the lamp has cooled down sufficiently (approx. 5 minutes) to be used again. If the problem persists, contact Technical Service.

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For further information please see the Operating Instruction of the Manufacturer.

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Data values are valid as output of the ø8mm optical fibre. With the ø5mm optical fibre, power and energy values increase by 46%.

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
The headrest can be more or less extracted from the backrest only by pulling or pressing it because it is automatically blocked by a clutch. The double joint for wide articulation can be adjusted by the knob (1):

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- In the case of the **standard** head (Fig. 12) rest turned in anti-clockwise sense, it will be unblocked; in clockwise sense it will be blocked;

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
- In case of optional version pneumatic "**Bios**" (Fig. 13); maintaining pressed the push button, the head rest is unblocked; leaving the push button the head rest is blocked.

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The right armrest can be lowered towards the front side chair to facilitate the patient's access. For doing it turn it forwards and downwards. For restoring it in the normal position make the opposite operation.

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Ordinary maintenance

The alimentation of the spray bottle

Before taking off the clean bottle, be sure that the LED indicator situated on the front of the dentist's instrument table has to be switched off, this guarantees that inside the bottle there is no pressure.

For taking off the bottle, turn it in anticlockwise sense.

To reinstall it repeat the inverse operation.

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WARNING: do not fill the bottle with saline solution. Fill the bottle only with water and, if necessary, disinfectant.

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For ensuring the safety of the device, it is recommended to replace the bottle every 5 years.

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If present, the Swipure filter must be replaced on a yearly basis. contact our technical support service.

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For keeping the safety of the device, it is recommended to execute, every one (1) year, the tests of conformity to the electrical safety.

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Cleaning and disinfection of the unit

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It is advised to clean up the external surfaces of the unit with hygienic liquid for a plastic material. The instruments tubes must be washed with water and neutral soap, then well dry. The upholstery of the patient chair and the stools have to be cleaned up with water and neutral soap.


Fig. 22 • Pct. 22 • Abb. 22

Bacinella (Opzionale se presente)

Cuspidor bowl (Optional, where provided)

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The cuspidor and the nozzles of the cup filler and washing can be removed for cleaning and disinfection pulling them. The cleaning of the small items must be done by a disinfectant liquid. Use vaseline oil to lubricate the O-rings in the nozzle connector.


Fig. 23 • Pct. 23 • Abb. 23

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The pedal battery must be recharged periodically. It can be charged by connecting the cable pedal to the patient chair (Fig. 43) with the pedal switch in the ON position. Charging time is approximately 4-5 hours.

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Product Identification

For any communication with the manufacturer or after sales centres, it must be always mentioned the equipment's model (mod.) and the serial number printed on the matriculation label situated on the arm of the installation.

Technical notes

The wiring diagrams and the documents for assistance and maintenance are given available in the qualified technical central services. The technical assistance on the equipment must be done exclusively by a high qualified technician.

Disposing of the device

Definition of RAEE

The symbol "REEE" means Refusal of Electrical and Electronically Equipment; it has to be considered as RAEE all electrical and electronically equipments for domestic or professional use that at the end of their life must be "scraped".

Disposing modality of the equipment

At the end of the usefull life, expected by the manufacturer in 10 years from the installation date, the equipment must be made unusable following this instructions:

- Untie the electrical, hydraulic and pneumatic connections of the equipment;
- Cut all electrical connections which go from the base of the equipment to different peripherals;
- Cut the cables of the secondary exit of the transformer as they cannot be restored.

When the equipment is unusefull, it has to be packed by recyclable material and as to be decided from the customer or the retailer, it has to be sent to private or public rubbish dumps equipped for the disposing of the REEE.

The customer or the retailer must take trace of the made "Disposing" of the equipment through the suitable structures, keeping the transport documents with precise motive's indication.

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Proper separate waste collection and compliance with the instructions above help prevent adverse effects on human health and the environment and promote the re-use and/or recycling of the materials with which the equipment is made. Unlawful disposal of the product results in the application of the administrative sanctions provided by the laws in force.

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Attention!! The symbol placed on the equipment, does not permit the disposing of the same with different ways from what we explained above.

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